

Gold Balances



GX/GF-A Series

1100 x 1mg
1620 x 1mg



FZ/FX-GD Series

320 x 1mg
520 x 1mg

Specifications	1mg Balances			
Model	GX-1003A GF-1003A	GX-1603A GF-1603A	FZ-300GD FX-300GD	FZ-500GD FX-500GD
Weighing Capacity	1100g	1620g	320g	520g
Resolution	0.001 g			
Calibration	GX-A Internal GF-A External		FZ-GD Internal FX-GD External	
Display Type	Clear reverse backlit LCD		Vacuum Fluorescent Display (VFD)	
Weighing Pan	128 mm x 128 mm		Ø130 mm	
Standard Interface	RS-232C and Quick USB		RS-232C	




GF-3000

3100g x 10mg



FZ/FX-GD Series

3200g x 10mg
5200g x 10mg



EK-610GD

600g x 10mg



EK-1202V

1200g x 10mg



EK-610V

610g x 10mg



HT-120

120g x 10mg



Specifications	10mg Balances						
Model	GF-3000	FZ-5000GD FX-5000GD	FZ-3000GD FX-3000GD	EK-610GD	EK-1202V	EK-610V	HT-120
Weighing Capacity	3100g	5200g	3200g	600g	1200g	610g	120g
Resolution	0.01g						
Calibration	External	FZ-GD - Internal FX-GD - External		External			
Display Type	Vacuum Fluorescent Display (VFD)			LCD With Backlit			
Weighing Pan	165mm x 165mm	ø150mm		ø110mm	127mm x 140mm	ø110mm	ø70mm
Standard Interface	RS-232C						N/A