

Hallmarking Balances



BM-5DH
Micro Analytical Balance
2.1g x 0.001mg / 5.2g x 0.01mg



GR-202
Semi-Micro Analytical Balance
42g x 0.01mg / 210g x 0.1mg



GF-3000
Precision Balance
3100g x 10mg

Touch Balances



HR-251AZ/A
Compact Analytical Balance
62g x 0.1mg / 252g x 1mg



FZ/FX-300GD
Precision Balance
3200g x 10mg

Specifications	Hallmarking Balances			Touch Balances	
Model	BM-5DH	GR-202	GF-3000	HR-251AZ/A	FZ/FX-300GD
Weighing Capacity	2.1g / 5.2g*	42g / 210g	3100g	62g / 252g	3200g*
Resolution	0.001mg / 0.01mg	0.01mg / 0.1mg	10mg	0.1mg / 1mg	10mg
Calibration	Automatic Self Calibration		External	HR-AZ Internal HR-A External	FZ-GD Internal FX-GD External
Display Type	Clear reverse backlit LCD display	LCD	Vacuum Fluorescent Display (VFD)	Clear Reverse Backlit LCD	Vacuum Fluorescent Display (VFD)
Weighing Pan Size	ø25mm	ø85mm	165mm x 165mm	ø90mm	ø150mm
Standard Interface	RS-232C and Quick USB		RS-232C		

*Smart range function: BM-5DH - The resolution will switch to 0.01 mg automatically when the display value exceeds 2.1g but will return to 0.001mg by pressing the RE-ZERO key.

HR-251AZ/A - The resolution will switch to 1 mg automatically when the display value exceeds 62 g but returns to 0.1 mg by pressing the RE-ZERO key.