

SS-Ki Series



FEATURES

- ▶ Rechargeable battery
- ▶ Detachable display unit
- ▶ 22 mm height LCD with backlight
- ▶ Large weighing pan
- ▶ Gross, Net, Tare function
- ▶ Comparator function
- ▶ Accumulative function
- ▶ Animal weighing

Options & Accessories:

Printers Remote Display

Specifications	SS-10Ki	SS-15Ki	SS-30Ki	SS-40Ki
Capacity	10 kg	15 kg	30 kg	40 kg
Readability	0.5 g	0.5 g	1 g	1 g
Repeatability (Standard deviation)	0.5 g	0.5 g	1 g	1 g
Linearity	± 1.5 g	± 1.5 g	± 3 g	± 3 g
Display	6 digit LCD with backlight (Character height 22 mm)			
Operating temperature/ humidity range	+ 10 °C to + 40 °C Less than 85% R.H. (no condensation)			
Power supply	AC adapter (100-240VAC 12VDC/0.5A) & Standard Ni-MH battery (7.2V/2Ah)			
Battery operating hours (on a full charge)	Approximately 40 hours (with backlight off)			
Interface	Standard RS-232C			
Weighing pan size	320 mm X 360 mm			
Net weight	Approximately 9 kg			



A&D Instruments India Private Limited

D-48, Udyog Vihar, Phase -V,
Gurugram- 122 016, Haryana
Tel: +91 -124-4715555, Email: info@aanddindia.in,
Web: www.aanddindia.in

Branch Offices:

Kolkata | Mumbai | Hyderabad | Chennai | Bengaluru

Dealer's Address