

PR-2000

Dot Matrix Printer



Features:

- Password function
- Editable header & footer
- Sample name & batch number can be inputted & printed
- Compliance with GLP/FDA/GMP Certification
- Print mode Direct, Weighing, Multi Weighing, (Tare, Gross, Net), (Multi T,G,N), Statistic, Sum, Mean, Formulation, Loss Weight

- RTC (Real-Time Clock) for printing date & time
- High speed dot matrix printing
- Easy setting and operation with LCD display
- Interval printing mode for printing weighing data at set intervals
- Compact & lightweight for portability & installation anywhere

Print Samples:

```

A&D
Model: GR-200
S/N: 14263411
ID: 0C/0020
Date: 2024/12/18
Start: 19:33:43
-----<STAT>-----
001 100.0007 g
002 100.0009 g
003 100.0009 g
004 100.0009 g
005 100.0007 g
006 100.0008 g
-----<RESULT>-----
n 6
Total 600.0048 g
max 100.0009 g
min 100.0007 g
diff 0.0002 g
s 0.00000 g
srel 0.00009 %
End: 19:34:51
Signature:
  
```

Statistic

```

A&D
Model: GR-200
S/N: 14263411
ID: 0C/0020
Date: 2024/12/18
Time: 19:24:42
-----
Tare: 49.9991 g
Gross: 100.0008 g
Net: 50.0017 g
-----
End: 19:25:03
Signature:
  
```

GTN

```

A&D
Model: GR-200
S/N: 14263411
ID: 0C/0020
Date: 2024/12/18
Time: 19:29:03
-----
No. Time T(g) G(g) N(g)
001 19:29:23 19.9997 49.9991 29.9994
002 19:30:02 49.9992 100.0003 50.0016
003 19:30:23 100.0007 200.0021 100.0014
Total 169.9996 350.0020 180.0024
-----
End: 19:30:32
Signature:
  
```

Multiple GTN

```

A&D
Model: FZ-500GD
S/N: 14263411
ID: 0C/0020
Date: 21/Oct/2024
Start: 11:56:21
-----<FORMULATION>-----
COMP001 -0.009 g
COMP002 -0.008 g
COMP003 -0.013 g
COMP004 -0.007 g
COMP005 0.000 g
Net: -0.037 g
-----
End: 11:56:41
Signature:
  
```

Formulation

```

A&D
Model: FZ-500GD
S/N: 14263411
ID: 0C/0020
Date: 21/Oct/2024
Start: 11:59:52
-----<LOSS WEIGHT>-----
Sample: 128.40 g
Lost: 9.67 g
Final: 118.73 g
-----
End: 12:00:12
Signature:
  
```

Loss Weight

SPECIFICATIONS	PR-2000
Printing method	Impact type dot matrix
Printing speed	Approx. 1.7 lines per second
Character size	1.7 mm (W) X 2.6 mm (H)
No. of characters printed	24 characters per line
Paper type	Normal (non-thermal)
Paper size	57 mm (W) X 50 mm (max diameter)
Operating temperature range	5°C to 45°C
Operating humidity range	10% to 80%RH
Interface	RS-232C D-Sub9 male
Power Supply	AC adapter - 12V (1.5A)
External dimensions	115 (W) x 190 (D) x 90 (H) mm
Weight	Approx. 600 g
Standard accessories	Ink ribbon (1 pc), plain printer paper (1 roll), AC adapter, RS-232C cable