

Wheelchair Weighing Scale

Easily accessible by people using wheelchairs

ADFSA-500KW-9090

FEATURES

- ▶ Tare function
- ▶ Large easy to read display
- ▶ Approved by Legal Metrology
- ▶ Inbuilt rechargeable battery
- ▶ Anti skid scale structure

**Colour options: Grey standard & black on order only.*



Front stopper



Side guard ensure safe loading



Large display



Low profile scale allows easy access



Specifications	ADFSA-500KW-9090
Weighing capacity	500 kg
Readability (d)	0.05 kg
MOC - Indicator	ABS Plastic
MOC - Structure	MS powder coated
Display	LCD with white backlight (character height 52 mm)
Operating Temperature / humidity range	+ 10 °C to + 40 °C Less than 85% R.H. (non-condensing)
Power supply	AC adapter (100-240V AC / 12V -500 mA) Rechargeable lead acid battery (6V/4Ah)
Battery operating hours (on a full charge)	Approximately 50 hours with backlight off
Platform size	900 mm X 900 mm
Net weight	Approximately 70 kg