


A Complete Weighing Solution for Jewellery Industry from **A&D Company, Ltd. Japan**



Hallmarking Balance



Gold Balance



Carat Balance



Silver Balance



Hallmarking Balances



BM-5DH

Micro Analytical Balance
2.1g x 0.001mg / 5.2g x 0.01mg



GR-202

Semi-Micro Analytical Balance
42g/210g x 0.01mg/0.1mg



GX/GF-3002A

Precision Balance
3200g x 10mg

Touch Balances



HR-251AZ/A

Standard-Level Analytical Balance
62g/0.1mg / 252g x 1mg



FZ/FX-3000GD


Precision Balance
3200g x 10mg

Specifications	Hallmarking Balances			Touch Balances	
Model	BM-5DH	GR-202	GX-3002A GF-3002A	HR-251AZ/A	FZ/FX-3000GD
Weighing Capacity	2.1g / 5.2g*	42g / 210g	3200g	62g / 252g	3200g
Readability	0.001mg / 0.01mg	0.01mg / 0.1mg	10mg	0.1mg / 1mg	10mg
Calibration	Automatic Self Calibration		GX-A Internal GF-A External	HR-AZ Internal HR-A External	FZ-Internal FX-External
Display Type	Reverse LCD with Backlit	LCD	Reverse LCD with Backlit		Vacuum Fluorescent Display (VFD)
Weighing Pan	ø25mm	ø85mm	165mm x 165mm	ø90mm	ø150mm
Standard Interface	RS-232C and Quick USB		RS-232C and Quick USB	RS-232C	


***Smart range function:** BM-5DH - The readability will automatically change to 0.01mg once the weight exceeds 2.1g. If needed, you can return the readability to 0.001mg by simply pressing the "RE-ZERO" key.

HR-251AZ/A - The readability will switch to 1 mg automatically when the display value exceeds 62 g, and will return to 0.1 mg by pressing the RE-Zero key.

Gold Balances



FZ/FX-GD Series
122g x 1mg (Only for FX-GD)
320g x 1mg
520g x 1mg




GX/GF-A Series
1100g x 1mg
1620g x 1mg


Specifications	1mg Balances				
Model	FX-120GD	FZ-300GD FX-300GD	FZ-500GD FX-500GD	GX-1003A GF-1003A	GX-1603A GF-1603A
Weighing Capacity	122g	320g	520g	1100g	1620g
Readability	1mg				
Calibration	External	FZ-GD Internal FX-GD External		GX-A Internal GF-A External	
Display Type	Vacuum Fluorescent Display (VFD)			Reverse LCD with backlit	
Weighing Pan	Ø130 mm			128mm x 128mm	
Standard Interface	RS-232C			RS-232C and Quick USB	



HT-120
120g x 10mg



EK-610GD
600g x 10mg




EK-610V
610g x 10mg




EK-1202V
1200g x 10mg



FZ/FX-GD Series
1220g x 10mg (only for FX)
3200g x 10mg
5200g x 10mg



GX/GF-3002A
3200g x 10mg



GX/GF-10202M
10200g x 10mg

Specifications	10mg Balances								
Model	HT-120	EK-610GD	EK-610V	EK-1202V	FX-1200GD	FZ-3000GD FX-3000GD	FZ-5000GD FX-5000GD	GX-3002A GF-3002A	GX-10202M GF-10202M
Weighing Capacity	120g	600g	610g	1200g	1220g	3200g	5200g	3200g	10200g
Readability	10mg								
Calibration	External					FZ-GD - Internal FX-GD - External		GX-A Internal GF-A External	GX-M Internal GF-M External
Display Type	LCD with Backlit				Vacuum Fluorescent Display (VFD)			Reverse LCD with Backlit	
Weighing Pan	ø70mm	ø110mm		127mmx140mm	ø150mm			165mmx165mm	270mmx210mm
Standard Interface	N/A	RS-232C						RS-232C and Quick USB	

www.aanddindia.in

Carat Balances




HR-AZ/A Series
310ct - 1260ct x 0.001ct


Specifications	Carat Balances		
Model	HR-250AZ/A	HR-150AZ/A	HR-251AZ/A
Weighing Capacity	1260ct / 252g	760ct / 152g	310ct / 1260ct 62g / 252g*
Readability	0.001ct / 0.1mg		0.001ct / 0.01ct 0.1mg / 1mg
Calibration	HR-AZ Internal HR-A External		
Display Type	Reverse LCD with backlit		
Weighing Pan Size	ø90mm		
Standard Interface	RS-232C		

***Smart range function:** The readability will switch to 0.01ct/ 1 mg automatically when the display value exceeds 310ct /62 g, and will return to 0.001ct/0.1 mg by pressing the RE-Zero/tare key.


Silver Balances




EK-6100V
6100g x 100mg




SS-KV Series
15kg x 200mg /
33kg x 500mg



ES-KV Series
15kg X 100mg /
36kg x 200mg



GX-M Series
32kg x 100mg



GX-L Series
62kg x 100mg

Specifications	Silver Balances						
Model	EK-6100V	SS-15KV	SS-30KV	ES-15KV	ES-30KV	GX- 32001M	GX- 62001L
Weighing Capacity	6100g	15kg	33kg	15kg	36kg	32kg	62kg
Readability	100mg	200mg	500mg	100mg	200mg	100mg	
Calibration	External					Internal	
Display Type	LCD With Back lit	Large LCD (22 mm Character Height)		Large LCD (24 mm Character Height)		Reverse LCD with Backlit	
Weighing Pan Size	127mm x 140mm	320mm x 360mm		300mm x 210mm		270mm x 210mm	384mm x 344 mm
Standard Interface	RS-232C					RS-232C & Quick USB	

Accessories



ED-1006



PR-2000

Options for EK-V Series EJ-13 Density Determination Kit



Free Extended Warranty

AMC & Calibration Services



CC-3722

NABL Accredited



Download Brochure

Dealer's Stamp

A&D Instruments India Private Limited
D-48, Udyog Vihar, Phase-V, Gurgaon-122016,
Haryana, India Ph: +91 124-4715555
Email: info@aanddindia.in Web: www.aanddindia.in

A&D
Discover Precision

A&D Company, Limited, Japan
3-23-14 Higashi-Ikebukuro, Toshima-Ku,
Tokyo 170-0013 Japan
Ph: +81 3-5391-6132 Web: www.aandd.jp

Branch Offices : Kolkata • Mumbai • Hyderabad • Chennai • Bangalore